



PO Box 1289 - 16 Simulaids Drive  
Saugerties NY 12477

Toll Free: (800)-431-4310 Fax: (845)-679-8996

[www.simulaids.com](http://www.simulaids.com)

[info@simulaids.com](mailto:info@simulaids.com)



**Child & Infant  
Airway Management Trainer  
Instructions for use**

**Thank you for buying an American made product.**

## **Recommendations:**

Child Accepts up to 4.0 mm. ET Tube oral and nasal.

Infant Accepts up to 2.5 mm. ET tube oral and nasal.

E.T. tube not included. NEW! Airway Bifurcation.

## **Setup:**

Your airway management trainer is ready to use when delivered.

**NO SETUP IS REQUIRED.**

## **Lubrication:**

Use water based lubricants only.

A tube of KY Jelly is included with the purchase of your airway management trainer. Place a coating of KY Jelly on the E.T. tube before inserting it into the trainer. Use of a petroleum-based lubricant will void your Simulaids warranty.

## **Cleaning:**

All Simulaids training devices can be cleaned with mild soapy water. Rinse thoroughly with clean water, let dry before storing.

## **Repair:**

Simulaids airway trainers can only be repaired by our factory. Any attempt to repair the trainer by the user will void our warranty.

## **WARRANTY**

Simulaids warrants this product to be free from any defect in materials and/or workmanship for a period of three years from the date of purchase, as evidenced by the date of invoice when the product was shipped to the end user. This warranty expressly does not cover abuse, accidental or purposeful damage, or any form of modification to the product. Simulaids reserves the right to either repair or replace affected parts or the entire unit, at their sole discretion, after investigating and reviewing the actual product and the damage. In most instances, a digital photo of the product in question showing the damage will help qualify a product for return to the factory. At no time will any product be accepted at the plant without proper return authorization issued by Simulaids.

Freight and Shipping charges are the sole responsibility of the end user. No product will be received with shipping charges due. Any product considered for warranty work must be identified by serial number and invoice number from the agency through whom the product was purchased. Without this information the product will not receive a return authorization number as required.

### **Return Goods Policy:**

Should it be necessary to return an item for any reason, contact our Customer Service Department to obtain a RETURN AUTHORIZATION NUMBER. Merchandise cannot be accepted without prior return authorization. Please refer to your invoice number when phoning in your request for returning merchandise. If needed, a shipping carton and packing instructions will be sent to you to facilitate return of the item.

### **Parts Numbers:**

#101-115 - Infant Trainer

Sh. Size 15" x 10" x 8" Sh. wt. 4 lb.


#101-125 - Child Trainer

Sh. Size 15" x 10" x 8" Sh. wt. 5 lb.

#PP0363 - KY Lubricant (4 oz.)

For 3 Pack Replacement Lungs/ Stomach for both the INFANT and CHILD contact customer service.

# Simulaids



PO BOX 1289 - 16 Simulaids Drive  
Saugerties NY 12477  
Toll Free: (800)-431-4310 Fax: (845)-679-8996  
[www.simulaids.com](http://www.simulaids.com)